




Portable GOSonic Toothbrush

Where Style Meets Function

Portable Toothbrush Instructions



General description

- A. Brush head
- B. Power button
- C. Battery cover
- D. Protective cover



Features

- Toothbrush vibration rate: Up to 30,000 strokes per minute.
- 2 Speeds.
- Replaceable toothbrush head.
- Battery use: One AAA type battery.



How to insert the battery

- Turn the battery cover counter-clockwise as ① and pull down the cover to remove it.
- Insert battery with polarity as directed.
- Replace battery cover and turn clockwise as ② until the battery cover is securely in place.

Batteries should be recycled or disposed of as per state and local guidelines.

How to change the replacement brush head

Turn brush head as ① and pull up as ②. Align the new brush head as ① and turn it as ②. Notch on neck should align with inner notches of toothbrush head to lock in place.



Caution

- Please stop using this product and contact your dentist if excessive bleeding occurs after use or bleeding continues to occur after one week of use.
- When brushing your teeth, please use the vibration force from the toothbrush itself. Do not force the bristles hard against the teeth or gums.
- Avoid high heat and direct sunlight.
- Do not use wet hand to replace the battery. It will cause the battery to malfunction.
- If the power is low, the vibration frequency of the product will be reduced or stopped. This is not a malfunction, please change the battery.
- Clean and dry brush and cover before storing.
- Store in a place out of reach of children.

GOSonic toothbrush one year limited warranty

Pop Dental will replace your originally purchased item, within one year of the shipping confirmation date, if it is found to be defective in materials or function. This warranty is valid for the original purchaser. Proof of purchase is required - even if your product is a gift. After 30 days from original purchase date, Pop Dental will require that the defective unit be returned. The customer is required to pay the postage to return the defective unit. Electronic units will be evaluated by our service center. If your unit is repairable it will be serviced and shipped back to you. If your unit is not repairable Pop Dental will cover the cost to reship and replace your unit with a new unit. If a unit is returned and is functioning it will be returned to you.

*This warranty gives you specific legal rights and you may also have other rights which vary from state to state.

How to obtain warranty service

Go to www.popdental.com/register. Once we receive your request, we will give you detailed instructions on how to process the replacement of your product.

Call US: 1 (631) 619-7122, Mon-Fri 8am-5pm, MST, or email us at customerservice@popdental.com

Warranty Exchanges for 3rd Party purchases (Items Not Purchased Directly Through Pop Dental): Please provide proof of purchase with a date from your third party vendor. Any warranty exchange will require an approval from Pop Dental.

What is not covered under warranty

Broken brush head tips are considered damage to the product and are not covered under warranty.

Damage caused by misuse - submersion in water, abuse, neglect, alterations, or dropping the product.

Normal wear and tear, including chips, scratches, abrasions, discoloration or fading.

Limitation of remedies

In no event shall Pop Dental be liable for any special, incidental or consequential damages based upon breach of warranty, breach of contract, negligence, tort, or any other legal theory. Such damages include, without limitation, loss of savings or revenue; loss of profit; loss of use; the claims of third parties including, without limitation, dentists and dental hygienists; and cost of any substitute equipment or services.